



AA-F0017

TULLGREN FUNNEL SET

For the study of ecology of soil by extracting motile organisms using the effects of heat and light. Also can be use for both dry and wet extraction of soil organism.

This apparatus comprises;

- 1 pc cylindrical lamp house 25cm L x 11cm dia with a BC lamp holder mounted on a.s.s support which can be moved up or down to control the surface temperature of the soil. Fitted with 1 m 3-core cable. The lamp house is bolted to a retort boss for use on a vertical support rod 10mm to 13mm dia
- 1 pc Lamp, BC 40W 240V a.c
- 1 pc plastic funnel 110mm dia, polySOL with fine netting, fitted with a 75mm L transparent PVC tubing
- 1 pc Plastic Funnel 110 mm dia polySOL with coarse netting, fitted with a 75mm L transparent PVC tubing
- 1 pc retort ring, metal approx 9cm int dia
- 1 pc beaker, pp, 100ml polySOL
- **N.B Stand (Base and Rod) not included**